

Laboratory Ultra Low Deep-Freezer, chest

TC 904-85



External Dimension:	W = 1560 mm D = 695 mm (plus handle and hinge) H = 890 mm	
Inside Dimension:	Utility Space 1 W = 1190 mm D = 440 mm H = 635 mm	Utility Space 2 W = 210 mm D = 440 mm H = 385 mm
Capacity:	368 I	
Temperature Range:	-60°C to -85°C	

(III. similar)

Housing

galvanized sheet steel with white coating. Cover double-walled with door lock, adhesive band. With 4 castors.

Interior space

special aluminium with round edges

Isolation

Cyclopentane insulation, without space, free of CFCS

Interior fittings Optional: white coated wire baskets



(III. similar)

Refrigerating unit

fully hermetically sealed fitted on vibration-absorbing mounts (ambient temperature 25°C), air cooled, low noise (55 dB(A)), energy saving compressor with high quality vaporisation system

High quality electronic temperature controller.

Temperature irregularities are viewed acoustically Battery backup for temperature display and alarm messages

Actual and set value display digital. Actual value permanently readable Set value digital adjustable by switch. Working range from -60°C to -85°C

tritec ® Gesellschaft für Labortechnik und Umweltsimulation mbH Hüttenstraße 9 D-30165 Hannover

Homepage
E-mailwww.tritec-klima.com
info@tritec-klima.dePhone++49-511/3523508Fax++49-511/3521715technischeÄnderungen vorbehalten





Laboratory Ultra Low Deep-Freezer, chest

TC 904-85

Defrost

manually

Electrical dates

Power supply230 V/50 Hz / single-phaseOptional: 230 V/ 60 HzPower consumption:400 WEnergy consumption10,2 KW (24 Std.)Power cable2,0 m with schuko plug

Packing details (palletized)

Dimensions:	approx. 160x80x120 cm
Net weight:	88 kg
Gross weight:	110 kg
Country of Origin:	European Union
Customs clearance code:	8418 3080

tritec ® Gesellschaft für Labortechnik und Umweltsimulation mbH Hüttenstraße 9 D-30165 Hannover

Homepagewww.tritec-klima.comE-mailinfo@tritec-klima.dePhone++49-511/3523508Fax++49-511/3521715technischeÄnderungen vorbehalten



Laboratory Ultra Low Deep-Freezer, chest

TC 904-85



tritec ... and life gets guaranteed!

Qualifications

DQ (Design Qualification)

Definition: Documented proof that the quality-related, GMP-related requirements has been adequately addressed in the design of equipment, including buildings, premises and auxiliary equipment

The user-requirement profiles (specifications) are documented and confirmed by us. On request, a specification can be created by us.

IQ (Installation Qualification)

Definition: Documented proof that critical equipment and systems have been delivered and installed in accordance with the set requirements and government regulations.

The IQ documentation is worked out by us especially for the delivered machine and is made available to you. The IQ documentation has to be carried out by the customer itself.

OQ (Operational Qualification)

Definition: Documented proof that critical equipment and systems in accordance with the set requirements in the whole operating range are working as intended in accordance with predetermined limits.

The OQ documentation is worked out by us especially for the delivered machine and is made available to you. The OQ documentation has to be carried out by the customer itself.

CQ (Calibration Qualification) according to DIN 13277:2022-05

Definition: Documented proof that critical measuring equipment in the intended range in accordance with predetermined tolerances operate reliably under current operating conditions

Verifying the temperature in the unloaded cooling unit (after reaching the steady state)

1 temperature on 3 measuring levels with 5 measuring points each

(Measurement with calibrated PT 1000 sensors). Test time 4 hours, then open door for 60 seconds. During this time, the limit values specified in DIN 13277:2022-05 must not be exceeded. Repeat the door opening after one hour.

The temperature measurements are carried out on our premises. The evaluation of the measurements, including graphical representation, is made in written form. The values must not exceed the limit values specified in DIN 13277:2022-05. (Other measuring methods possible on request)

PQ (Performance-Qualification) according to DIN 13277:2022-05

Definition: Documented proof that critical equipment and systems in accordance with the set requirements in the whole workspace under current working conditions (with product) provide the requested services. The calibration described above is carried out under real conditions on site. Optionally, the measurement can be can be carried out in a loaded or unloaded state. The measurement evaluation, including graphical representation, is made in written form. The values must not exceed the limits specified in DIN 13277:2022-05. exceeded. (Other measuring methods possible on request)

Homepage
E-mailwww.tritec-klima.com
info@tritec-klima.dePhone++49-511/3523508Fax++49-511/3521715technischeÄnderungen vorbehalten

